Product:

Radix® fiber post

C0613 – C0614

Ref doc.:

F19 03 40.X/ 04 / 2005 updated 11/2009

PAGE

GB	DIRECTIONS FOR USE (Radix® fiber post)	2
F	PROTOCOLE D'UTILISATION (Radix® fiber post)	6
D	GEBRAUCHSANLEITUNG (Radix® fiber post)	11
I	MODO DI USO (Radix® fiber post)	16
E	MODO DE UTILIZACION (Radix® fiber post)	21
NL	GEBRUIKSVORSCHRIFT (Radix® fiber post)	26
DK	BRUGSANVISNING (Radix® fiber post)	31
FIN	KÄYTTÖOHJEET (Radix® fiber post)	36
S	BRUKSANVISNING (Radix® fiber post)	41
P	PROTOCOLO DE UTILIZAÇÃO (Radix® fiber post)	46
GR	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ ΓΙΑ (Radix® fiber post)	51
POL	INSTRUKCJA OBSŁUGI (Radix® fiber post)	56

For Dental Use Only**DIRECTIONS FOR USE (Radix® fiber post)****0) COMPOSITION**

Quartz fibers 60% volume
Epoxy resin 40% volume

1) INDICATIONS FOR USE

These instruments are to be used only in a clinical or hospital environment, by qualified users.

Application field

In case of insufficient residual tooth substance (< 4 mm), the Radix fiber is needed to support the coronal restoration.

2) CONTRAINDICATIONS

Insufficient residual dentin: at least 2 mm of tooth structure is required around the preparation.

3) WARNINGS

None known.

4) PRECAUTIONS

The post must be cleaned with alcohol before insertion in the canal.

Avoid reprocessing the post more than once because numerous reprocessing cycles might reduce the useful life of the device and damage the marking.

Avoid touching the posts with your fingers after cleaning.

Shortening the post should be done outside of the mouth.

Dentsply Maillefer recommends the use of a rubber dam.

5) ADVERSE REACTIONS

In the present technical state, no adverse reaction has been reported so far.

6) STEP-BY-STEP INSTRUCTIONS

1. Root canal preparation. Ensure that sufficient space exists to place the post. Remove residual gutta percha out of the root canal.
2. Select correct Radix fiber size, according to anatomical situation of tooth using the radiograph and the information provided below. Select the Largo® Peeso Reamer and Radix fiber Precision Drill corresponding to the selected Radix fiber size.

Radix fiber

	No 1	No 2	No 3	No 4	No 5	No 6
A.	0.8	0.8	1.0	1.0	1.3	1.3
B.	1.35	1.47	1.67	1.83	2.04	2.22
C.	No 1	No 1	No 2	No 2	No 3	No 3
D.	No 1	No 2	No 3	No 4	No 5	No 6

A. Tip diameter [mm]

B. Head diameter [mm]

C. Largo® Peeso Reamer

D. Radix fiber Precision Drill

3. Determine the length of the coronal build-up which will account for $\frac{1}{3}$ of the final post length. Remove the root canal obturation with the selected Largo® Peeso Reamer (rotation speed 800 – 1200 rpm) as deep as necessary for the post to be inserted at $\frac{2}{3}$ of its length. At least 4 mm of root canal filling should remain in the apical region. In case of curved canals, the length has to be reduced.
4. Drill the canal with the selected Radix fiber Precision Drill (rotation speed 1000 – 1200 rpm).
5. Check that the post is properly seated in the canal.
6. Shorten the post to its final length with a diamond disc outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.
7. Clean the post with alcohol.
8. Apply etchant in the canal and to the exposed dentin for 15 seconds. (e.g. Dentsply DeTrey® Conditioner 36%; *). Rinse for 10 sec. Dry gently with paper points but leave surface moist. Do not put the etchant in contact with the gingival. If it is not possible, prefer the solution of a self-etching primer.
9. Mix the primer and its self-cure activator (e.g. Dentsply Prime&Bond® NT™ Dual Cure; *). Apply 2 coats of the preparation with a brush (e.g. Dentsply DeTrey® Applicator Tips) in the root canal and leave for 20 seconds. Remove excess using paper points and gently air dry all surfaces. Enamel

and dentin surfaces should have a uniform, glossy appearance, if not, repeat the application.

10. Apply a single coat of primer to the post. Gently air dry for 5 seconds and light-cure the post for 10 to 20 seconds outside the mouth.
11. Mix the cement(e.g. Dentsply Calibra™; *) and apply it in the canal and optionally to the post using a Lentulo spiral filler (e.g. Dentsply Maillefer, REF A 0022) or a syringe tip. Seat the post immediately. Remove excess with appropriate instruments. Light-cure for 40 to 60 seconds, applying gentle pressure on the post with the tip of the curing light optic.
12. Apply 2 coats of bonding agent (e.g. Dentsply Prime&Bond® NT™ Dual Cure; *) to the exposed post, cement and coronal areas and leave for 20 seconds. Remove excess with a jet of air. Light-cure for 30 seconds.
13. Directly model the build-up using a composite core build-up material (e.g. Dentsply Core X™; *)

7) DISINFECTION, CLEANING AND STERILIZATION

Reprocessing procedure for dental instruments and implantable radicular devices

Foreword

For hygiene and sanitary safety purposes, all instruments not marked "sterile" must be cleaned, disinfected and sterilized before each usage to prevent any contamination. This concerns the first use as well as the subsequent ones.

Area of application

Disinfection and sterilisation before first usage and reprocessing procedures concerning:

A1. Instruments

Cutting instruments, (hand and engine driven) such as:

- Endodontic instruments (files, broaches, reamers, enlargers, endodontic burs, ultrasonic inserts);
- Rotary cutting instruments (Diamond burs, tungsten carbide burs, stainless steel drills, carbon steel burs).

Root canal filling instruments (Pluggers, spreaders, compactors).

Supports, kits and instrument organisers.

Hand instruments and clamps.

A2. Implantable devices

Dentinal and radicular posts made of steel, titanium and glass fibers.

Supports, kits and organiser systems for posts.

A3. Contra angle

B. Filling material: Only chemical disinfection (no sterilisation)

Gutta percha, Thermafil obturation devices.

Exclusion

- Equipment such as Motors, Apex locators and other devices with reprocessing procedures included in the individual Direction for Use.
- MTA, Glyde, TopSeal.

General recommendation

- Use only a disinfecting solution which is approved for its efficacy (VAH/DGHM-listing, CE marking, FDA approval) and in accordance with the DFU of the disinfecting solution manufacturer. For all metal instruments, it is recommended to use anticorrosion disinfecting and cleaning agents.
- For your own safety, please wear personal protective equipment (gloves, glasses, mask).
- The user is responsible for the sterility of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments where applicable after sterility.
- Limitations and restrictions on reprocessing:
The individual DFU indicates if the useful life of a device might be reduced by the number of reprocessing cycles. Furthermore, the appearance of defects such as cracks, deformations (bent, twisted), corrosion, loss of colour coding or marking, are indications that the devices are not able to fulfil the intended use with the required safety level.
- Single use marked instruments are not approved for re-use.
- The water quality has to be convenient to the local regulations especially for the last rinsing step or with a washer-disinfector.
- Tungsten carbide burs, plastic supports, hand instruments and NiTi instruments are degraded by Hydrogen Peroxide (H_2O_2) solution.
- NiTi Instruments are degraded if immersed more than 5 minutes in a solution of NaOCl at more than 5%.
- Instruments made of aluminium are degraded in presence of caustic soda solutions with mercury salt. Do not use acid ($pH < 6$) or alkaline ($pH > 8$) solutions.
- The washer-disinfector is not recommended for instruments made of aluminium, tungsten carbide or carbon steel.

Step-by-step procedure

A. Devices

		A3. Contra Angle	
		A2. Implantable devices	
		A1. Instruments	Following uses
			First use
1.	Disassembling	- Disassemble the device, if required	Silicone stops have to be removed.
2.	Pre-Disinfection	- Soak all instruments immediately after use in a detergent and disinfecting solution combined with proteolytic enzyme if possible.	<ul style="list-style-type: none"> - Follow instructions and observe concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or others defects on instruments). - The disinfecting solution should be aldehyde free (to avoid blood impurities fixation) and without di- or triethanolamines as corrosion inhibitor. - Do not use disinfecting solutions containing Phenol or any products which are not compatible with the instruments (See general recommendations). - For visible impurities observed on instruments a pre-cleaning is recommended by brushing them manually with soft material.
3.	Rinsing	- Abundant rinsing (at least 1 min)	<ul style="list-style-type: none"> - Use quality water in accordance with local regulations. - If a pre-disinfectant solution contains a corrosion inhibitor, it is recommended to rinse the instruments just before the cleaning.
4a.	Automated Cleaning with washer-disinfector	<ul style="list-style-type: none"> - Place the devices in a kit, support or container to avoid any contact between instruments or posts. - Put them in the washer-disinfector (Ao value > 3000 or, at least 5 min at 90 °C). 	<ul style="list-style-type: none"> - Discard any instruments with large obvious defects (broken, bent). - Avoid any contact between instruments or posts when placing in the washer-disinfector use kits, supports or container. - Follow instructions and observe concentrations given by the manufacturer (see also general recommendations). - Use only approved washer-disinfector according to EN ISO 15883, maintain and calibrate it regularly.
OR			
4b.	Manual Cleaning or assisted by an ultrasonic device	<ul style="list-style-type: none"> - Place the devices in a kit, support or container to avoid any contact between instruments. - Immerse in the disinfecting solution with cleaning properties, assisted by an ultrasonic device if suitable. 	<ul style="list-style-type: none"> - No visible impurities should be observed on the instruments. - Discard any instruments with large obvious defects (broken, bent, and twisted). - Follow instructions and observe concentrations and time given by the manufacturer (see also general recommendations). - The disinfecting solution should be aldehyde free and without di- or triethanolamines as corrosion inhibitor.
5.	Rinsing	- Abundant rinsing (at least 1 min)	<ul style="list-style-type: none"> - Use quality water in accordance with local regulations. - If a disinfecting solution contains a corrosion inhibitor, it is recommended to rinse the instruments just before the autoclaving. - Dry on a single use non-weaved cloth, or with a drying machine or filtered compressed air.
6.	Inspection	<ul style="list-style-type: none"> - Inspect devices and sort out those with defects. - Assemble the devices (stops) 	<ul style="list-style-type: none"> - Dirty instruments must be cleaned and disinfected again. - Discard instruments which show any deformations (bent, twisted), damages (broken, corroded) or defects (loss of colour coding or marking) affecting the resistance, the safety or the performance of the instrument or posts. - Protect carbon steel bur with corrosion inhibitor before packaging. - For Contra Angle: lubricate the device with an adequate spray before packaging
7.	Packaging	<ul style="list-style-type: none"> - Place the devices in a kit, support or container to avoid any contact between instruments or posts and pack the devices in "Sterilisation pouches". 	<ul style="list-style-type: none"> - Avoid any contact between instruments or posts during sterilization. Use kits, supports or containers. - Check the validity period of the pouch given by the manufacturer to determine the shelf life. - Use packaging which are resistant up to a temperature of 141°C (286°F) and in accordance with EN ISO 11607.
8.	Sterilization	- Steam sterilisation at: 134 °C / 273°F during 18 min.	<ul style="list-style-type: none"> - The instruments, posts and the plastic supports must be sterilized according to the packaging labelling. - Use only autoclaves that are matching the requirements of EN 13060, EN 285. - Use a validated sterilisation procedure according ISO 17665 - Respect the maintenance procedure of the autoclave device given by the manufacturer. - Use only this recommended sterilization procedure. - Control the efficiency (packaging integrity, no humidity, colour change of sterilisation indicators, physico-chemical integrators, digital records of cycles parameters). - Traceability of procedure records
9.	Storage	- Keep devices in sterilization packaging in a dry and clean environment	<ul style="list-style-type: none"> - Sterility cannot be guaranteed if packaging is open, damaged or wet. - Check the packaging and the medical devices before using them (packaging integrity, no humidity and validity period).

B. Filling material

<i>Operation</i>	<i>Operating mode</i>	<i>Warning</i>
1. Disinfection	- Immerse the obturation devices in NaOCl (2,5% at least) during 5 mn at ambient temperature.	- Do not use disinfecting solutions containing Phenol or any products which are not compatible with the treated filling material. (See general recommendation).

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